

# Kleargene XL

## Quick guide: Blood DNA extraction kit



### Cell lysis & DNA binding

- 30 mL buffer B1 ● + 10 mL whole blood, vortex
- Incubate for 5 min at room temperature, vortex
- Centrifuge (3,000 x g / 2 min), discard supernatant

### 1st buffer C1 wash

- + 30 mL buffer C1 ●, vortex
- Centrifuge (3,000 x g / 2 min), discard supernatant

### 2nd buffer C1 wash

- + 30 mL buffer C1 ●, vortex
- Centrifuge (3,000 x g / 2 min), discard supernatant

### 1st buffer A1 wash

- + 30 mL buffer A1 ●, vortex
- Centrifuge (3,000 x g / 2 min), discard supernatant

### 2nd buffer A1 wash

- + 30 mL buffer A1 ●, vortex
- Centrifuge (3,000 x g / 2 min), discard supernatant

### Buffer W1 wash

- + 30 mL buffer W1 ●, vortex
- Centrifuge (3,000 x g / 2 min), discard supernatant

### Ethanol wash

- + 30 mL 100% ethanol, vortex
- Centrifuge (3,000 x g / 2 min), discard supernatant

### Silica drying

- Invert open tube over tissue paper for 5 min
- Dry open tube in 55°C fan oven for 30 min

### DNA elution

- Add appropriate volume of pre-warmed buffer (@ 55°C) E1 ●, vortex
- Incubate tube in 55°C fan oven 15 min, vortex
- Centrifuge (3,000 x g / 5 min @ 4°C), aspirate off DNA solution

This protocol can be scaled to the appropriate volume of blood by reducing/increasing the volumes of wash buffers used in proportion; this protocol assumes that 10 mL of blood is to be extracted.

**Need help, please contact our technical support team on**  
 +44 (0)1992 470 757 or email us at [tech.support@lgcgenomics.com](mailto:tech.support@lgcgenomics.com)

[www.lgcgenomics.com](http://www.lgcgenomics.com)

Germany • United Kingdom • United States