

Kleargene

Quick guide: Kleargene plant 384-well DNA extraction kit



Sample disruption

- Homogenise samples, + 75 μ L buffer L1 ●, incubate (55°C / 10 min)
- + 150 μ L buffer C1 ●, mix, incubate (55°C / 10 min)
- Centrifuge (3,000 x g / 2 min) to clear lysate

DNA binding

- Place filter plate on top of collection plate / reservoir
- Transfer 100 μ L cleared lysate to filter plate wells, incubate at room temperature for 2 min
- Centrifuge (3,000 x g / 2 min) then clear collection plate

1st wash

- + 75 μ L buffer A1 ● to filter plate wells, incubate at room temperature for 2 min
- Centrifuge (3,000 x g / 2 min) then clear collection plate

2nd wash

- + 37.5 μ L buffer W1 ● to filter plate wells, incubate at room temperature for 2 min
- Centrifuge (3,000 x g / 2 min) then clear collection plate

Ethanol wash

- + 37.5 μ L 100% EtOH to filter plate wells, incubate at room temperature for 2 min
- Centrifuge (3,000 x g / 2 min) then clear collection plate

Drying & elution

- Incubate filter plate (55°C / 10 min), then place on an elution plate
- + 37.5 μ L buffer E1 ● (@ 55°C) to filter plate wells, incubate (55°C / 5 min)
- Centrifuge (3,000 x g / 2 min)

Need help, please contact our technical support team on
+44 (0)1992 470 757 or email us at tech.support@lgcgenomics.com

www.lgcgenomics.com

Germany • United Kingdom • United States