

# Kleargene

## Quick guide: Plant silica slurry DNA extraction kit



### Sample disruption

- Homogenise samples, + 150 µL buffer L1●, incubate (55°C / 30 min)
- Centrifuge (3000 x g / 2 min) to pellet debris
- Transfer as much cleared lysate as possible into a new plate

### DNA binding

- + 300 µL well mixed buffer B1●, mix, incubate for 5 min at room temperature
- Centrifuge (3000 x g / 2 min) to pellet silica, then decant supernatant

### 1st wash

- + 300 µL buffer C1●, mix, incubate for 5 min at room temperature
- Centrifuge (3000 x g / 2 min) then decant supernatant

### 2nd wash

- + 300 µL buffer A1●, mix, incubate for 5 min at room temperature
- Centrifuge (3000 x g / 2 min) then decant supernatant

### 3rd wash

- + 150 µL buffer W1●, mix, incubate for 5 min at room temperature
- Centrifuge (3000 x g / 2 min) then decant supernatant

### Ethanol wash

- + 150 µL EtOH, mix, incubate for 5 min at room temperature
- Centrifuge (3000 x g / 2 min) then decant supernatant

### Drying & elution

- Incubate (55°C / 30 min), then place on a new collection plate
- + 150 µL buffer E1● (@ 55°C), mix, incubate (55°C / 10 min)
- Centrifuge (3000 x g / 2 min) then carefully aspirate off DNA solution

**Need help, please contact our technical support team on**  
 +44 (0)1992 470 757 or email us at [tech.support@lgcgenomics.com](mailto:tech.support@lgcgenomics.com)

[www.lgcgenomics.com](http://www.lgcgenomics.com)

Germany • United Kingdom • United States