Efficient high-throughput SNP discovery and genotyping using normalised genotyping by sequencing (nGBS)

Abstract

Large-scale genotyping programmes require efficient tools for high-density SNP discovery. The commonly used technique involving whole genome shotgun sequencing of parental genotypes is often prohibitively expensive and resource intensive. Here we present an optimised, self-tuning <u>genotyping</u> by sequencing (GBS) method, called normalised GBS (nGBS), which efficiently reduces the genome complexity of any species to a few hundred thousand loci across the complete genome. No prior knowledge of the genome is required, and repetitive regions are avoided. The selected loci are reliably sequenced across large sample batches using <u>next generation sequencing</u>. Simultaneous SNP discovery and genotyping is carried out for between forty-eight and several thousand samples without prior assay design effort and bias.



Efficient high-throughput SNP discovery and genotyping using normalised genotyping by sequencing (nGBS)

Background

The majority of most recent projects in industry and academia involving large scale breeding programmes, trait mapping, germplasm characterisation and quality control now rely on the tools offered by molecular genetic markers (He et al. 2014). Traditionally, microsatellite markers have been used in such applications but, in the last decade, single nucleotide polymorphism (SNP) markers have become increasingly important. SNP markers are particularly efficient due to their specificity (as they target a single base pair difference) and to the availability of simple and robust genotyping methods.

Development of SNP markers is, however, relatively expensive and time-consuming, especially when very dense genetic maps are necessary. Although high-throughput sequencing costs are decreasing rapidly, the effort required for sequencing and data analysis is still prohibitive. This is a particular issue when performing whole genome shotgun (WGS) sequencing of large genomes, where many genotypes need to be analysed to fully cover the genetic diversity within the species.

Furthermore, the data analysis involved (de novo assembly, variant calling) in comparing the different varieties and developing accurate SNP PCR assays can be very time consuming and require expensive computational resources.

Traditional genotyping by sequencing

Instead of analysing the complete genome sequence of each sample, one can focus on reducing the complexity of the genome sequence to certain regions or features which can be cost-efficiently sequenced for a large

number of samples. One way to achieve this. with little or no prior knowledge of the genome in question, is to use restriction enzyme digestion of the genomic DNA. This involves the use of one or two restriction enzymes, and results in a reproducible set of fragments for each sample. These non-randomly selected DNA pieces represent genomic loci, and typically cover between several tens of thousands and a few hundred thousand loci in the genome. Sample fragments are then subjected to highly multiplexed sequencing, producing up to several million sequence reads for each sample on standard next generation sequencing instruments. The resulting sequence data is analysed in order to determine the polymorphic loci (typically every 10th to 1,000th of the covered loci) and subsequently the genotype of each polymorphic locus for each sample.

There are several protocols for this technique, with the following being the most common in the literature:

- RAD-Seq (Restriction Associated DNA sequencing) – single enzyme digestion combined with random fragmentation, sequencing adapter ligation (Baird et al. 2008)
- GBS single enzyme digestion, size selection and inline barcode addition with custom PCR protocol; typically using the ApeKI enzyme (Elshire et al. 2011)
- ddRAD-Seq double enzyme digestion, size selection and inline barcode addition with custom PCR protocol; often PstI-MspI (Poland et al. 2012)

Efficient high-throughput SNP discovery and genotyping using normalised genotyping by sequencing (nGBS)

The number of genomic loci covered per sample for a given method can be theoretically estimated if the average guanosine-cytosine (GC) content is known for the genome. However, there are other technical factors affecting the final restriction fragment profile, such as methylation pattern and fragment size selection during library preparation as well as clustering bias on the sequencing flow cell. Sequence diversity across the samples also influences the restriction fragment pattern; diversity is high in ecological studies of wild populations and low in inbred populations. This latter sequence diversity factor also influences the main outcome of the experiment, which is the number of reliably covered polymorphic loci across a majority of the samples.

The optimised normalised genotyping by sequencing method

When dealing with a new species, it is usually necessary to perform extensive experimental trials involving several enzyme combinations and various protocol modifications. The restriction enzyme(s) used are usually chosen to be sensitive to the most abundant types of methylation in the target genome and thus typically avoid repetitive regions (mainly retrotransposons). Some high copy number DNA in the sample (e.g. mitochondrial DNA, chloroplasts and ribosomal repeats) are not methylated, and can result in a reduced number of sequencing reads from the lower copy number regions.

In order to simplify the protocol development stage and to reduce the sequencing required due to the presence of non-methylated high copy number DNA, LGC Biosearch Technologies has developed a self-tuning protocol applicable to any kind of genome, regardless of size, methylation pattern, and type and abundance of repetitive elements. This method is called normalised genotyping by Sequencing (nGBS) and an overview is given in figure 1. The normalisation step is adapted from hybridisation kinetics resulting in reduction of abundant fragments. nGBS uses the MsII restriction enzyme to produce blunt end fragments and an enzyme treatment in a subsequent normalisation step. Although MsII digestion is inhibited by overlapping CNG and CHH methylation (common methylation types in plants), the number of fragments produced is generally too high for cost-efficient sequencing. Therefore, the normalisation step is important to reduce any remaining high copy number fragments.

Based on prior genome knowledge and the availability of GBS data for comparison, a suitable sequencing protocol is chosen. This can range from single-sided 75 bp sequencing (at minimal cost per sample) for species with a genome of low complexity and known sequence, to paired 150 bp sequencing which is the protocol of choice when working with species for which no reference sequence is available. Up to 384 samples are multiplexed and sequenced on an Illumina NextSeg flow cell; this offers the flexibility of generating a few hundred thousand reads to a few million read pairs per sample. Generally, an average depth of 1.5 million read pairs per sample is sufficient for a well-tuned protocol.

The appropriate data analysis is again dependent on prior knowledge of the genome and availability of previous data:

• With an available reference sequence or a well- saturated reference of GBS tags, the

Efficient high-throughput SNP discovery and genotyping using normalised genotyping by sequencing (nGBS)

newly acquired sequence data is aligned, observed alleles called, and genotypes determined for each sample according to the read counts seen for each allele.

 When no reference sequence, GBS tags or fragment cluster collection is available the newly acquired sequence data is clustered to generate a reference. For 75 bp protocols, this is performed using all unique sequence tags from the first 64 bases as input. For paired 150 bp protocols, where the typical insert sizes are below 270 bp, all unique sequences from paired end overlap consensus fragments are used. The polymorphic tag cluster are determined. Genotypes are called for each sample based on the read counts for each allele.

Either data analysis pathway results in a genotype table with the samples as columns, and the SNP locus as rows. Filters are then applied to remove low-quality loci, typically resulting in a set of 5,000-50,000 polymorphic loci with high-quality genotype calls in ≥75% of the samples (dependent on population diversity and sequence depth per sample). The entire data analysis process can be carried out with publicly available GBS analysis packages e.g., Stacks, TASSEL or PyRAD, or an Biosearch Technologies-developed pipeline optimised for nGBS data.

Low sampling rate in sequencing of the libraries or the presence of existing polymorphisms in the restriction sites (presence/absence variants) result in missing genotype data. As certain downstream applications (e.g. GWAS, QTL) typically expect complete datasets, missing data imputation (e.g. Swarts et al. 2014) can then be used to correct and complete the dataset, a process that is facilitated by previous knowledge of the organism (published tag sets) and the population studied.

Application of nGBS

The method proposed here has following applications: Co-dominant SNP genotyping, Genetic diversity analysis, Phylogenetic study, Genetic mapping, QTL mapping, Population genetics, Mapping to whole genome sequence, Discovery of (rare) SNP variants and Genome wide association mapping.

The technique nGBS have been applied to genotype range of complex species including large genomes- Sitka spruce - Picea sitchensis (20Gb) for genetic diversity analysis (Byrne, T.,et al. 2022). Phylogenetic studies of Macaranga species was performed using nGBS (Dixit, N. M.et al. 2023), Population genomics study for large genome as Pyrenean larkspur -Delphinium montanum (Salvado, P et al. 2022) was carried out using nGBS, QTL mapping of peanut -Arachis hypogaea L. to identify gene associated with resistance to a diseases (Han S, et al. 2018). The large genomes as Scots pine- Pinus silvestris (22Gb), and Tulip- Tulipa (34Gb), Wheatgrass - Thinopyrum intermedium (13Gb), Timothy grass - Phleum pratense (12Gb). were genotyped using nGBS.

Efficient high-throughput SNP discovery and genotyping using normalised genotyping by sequencing (nGBS)

Conclusions

Biosearch Technologies has established the improved nGBS protocol as a very efficient tool in SNP discovery for marker development. The identified polymorphic locus sequences can be used directly to develop PCR-based SNP assays. The method can also be used directly for routine high- throughput genotyping, as the cost per sample is lower than SNP microarray genotyping ('SNP chips') with an optimised nGBS protocol. Additionally, the experimental design is faster with a lower associated cost and there is no design bias, which is especially important when working in species with very high polymorphism rates e.g. maize (see Elshire et al. 2011).

All these features make the nGBS method an ideal solution for rapid and high-quality SNP discovery and genotyping in any organism, with a reasonable price tag.




Figure 1. The nGBS protocol

Efficient high-throughput SNP discovery and genotyping using normalised genotyping by sequencing (nGBS)

References

- Dixit NM, Zirpel M, Slik JWF, Jamsari J, Weising K and Guicking D (2023) Biogeography of the Sunda Shelf revisited: Insights from Macaranga section Pruinosae (Euphorbiaceae). Front. Ecol. Evol. 10:1049243. <u>https://doi.org/10.3389/ fevo.2022.1049243</u>
- He J, Zhao X, Laroche A, Lu Z-X, Liu H and Li Z (2014) Genotyping-by-sequencing (GBS), an ultimate marker-assisted selection (MAS) tool to accelerate plant breeding. Front. Plant Sci. 5:484. <u>http://dx.doi.</u> org/10.3389/fpls.2014.00484
- Baird NA, Etter PD, Atwood TS, Currey MC, Shiver AL, Lewis ZA, et al. (2008) Rapid SNP Discovery and Genetic Mapping Using Sequenced RAD Markers. PLoS ONE 3(10): e3376. <u>http://dx.doi.org/10.1371/</u> journal.pone.0003376
- Elshire RJ, Glaubitz JC, Sun Q, Poland JA, Kawamoto K, Buckler ES, et al. (2011) A Robust, Simple Genotyping-by-Sequencing (GBS) Approach for High Diversity Species. PLoS ONE 6(5): e19379. <u>http://dx.doi.org/10.1371/journal.pone.0019379</u>
- Poland JA, Brown PJ, Sorrells ME, Jannink J-L (2012) Development of High-Density Genetic Maps for Barley and Wheat Using a Novel Two-Enzyme Genotypingby-Sequencing Approach. PLoS ONE 7(2): e32253. <u>http://dx.doi.org/10.1371/journal.pone.0032253</u>
- Swarts, K., H. Li, J. A. Romero Navarro, D. An, M. C. Romay, S. Hearne, C. Acharya, J. C. Glaubitz, S. Mitchell, R. J. Elshire, E. S. Buckler, and P. J. Bradbury. 2014. Novel Methods to Optimize Genotypic Imputation for Low-Coverage, Next-Generation Sequence Data in Crop Plants.

f in @LGCBiosearch

biosearchtech.com

All trademarks and registered trademarks mentioned herein are the property of their respective owners. All other trademarks and registered trademarks are the property of LGC and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details. No part of this publication may be reproduced or transmitted in any form or by any means, electronic or mechanical, including photocopying, recording or any retrieval system, without the written permission of the copyright holder. © LGC Limited, 2023. All rights reserved. GEN/1133/LH/0723

For Research Use Only. Not for use in diagnostic procedures.

- Plant Genome 7. <u>http://dx.doi.org/10.3835/</u> plantgenome2014.05.0023
- Han S, Yuan M, Clevenger JP, Li C, Hagan A, Zhang X, Chen C and He G (2018) A SNP-Based Linkage Map Revealed QTLs for Resistance to Early and Late Leaf Spot Diseases in Peanut (Arachis hypogaea L.). Front. Plant Sci. 9:1012. doi: 10.3389/fpls.2018.01012
- Byrne, T., Farrelly, N., Kelleher, C., Hodkinson, T. R., Byrne, S. L., & Barth, S. (2022). Genetic Diversity and Structure of a Diverse Population of Picea sitchensis Using Genotyping-by-Sequencing. Forests, 13(9), 1511. https://doi.org/10.3390/f13091511
- Salvado, P., Aymerich Boixader, P., Parera, J., Vila Bonfill, A., Martin, M., Quélennec, C., ... & Bertrand, J. A. (2022). Little hope for the polyploid endemic Pyrenean Larkspur (Delphinium montanum): Evidences from population genomics and Ecological Niche Modeling. Ecology and Evolution, 12(3), e8711.

