



LGC

For Research Use Only. Not for use in diagnostic procedures.

All-inclusive KASP service guidance notes

### Contents

1.	Introduction	. 2				
2.	Project workflow	. 3				
3.	Sample and sequence submission	. 3				
4.	Data analysis and results	. 3				
5.	Retention and return of project samples	. 4				
6.	Useful links and contacts	. 4				
Арр	pendix A: Forms for services	. 5				
Арр	pendix B: DNA submission requirements	. 6				
Арр	Appendix C: Tissue submission requirements7					

### 1. Introduction

These guidance notes have been developed to assist with the submission of samples and supplementary information to LGC, Biosearch Technologies<sup>™</sup> for an all-inclusive services project.

The all-inclusive DNA extraction and KASP<sup>™</sup> genotyping service offering includes KASP genotyping assay design, DNA extraction and KASP genotyping services, and is delivered from our global state-of-the-art laboratories.

The complete service package from Biosearch Technologies includes:

- BioArk<sup>™</sup> Sample Collection Kit (for plant and aquaculture samples)
- DNA extraction service
- Custom KASP assay design
- KASP genotyping service and support

The complete service package is available for the following industries\*:

- Livestock (Animal samples (including buccal swabs, blood, ear punch, please enquire for other sample types) and DNA)
- Aquaculture (Animal samples and DNA)
- Plant (Plant samples or extracted DNA only)

\*Applicable surcharges may be incurred dependent on the sample type.

Please note that our all-inclusive service is based on *in-silico* designed KASP genotyping assays. The service does not include wet lab validation of assays or re-design as standard (options are available, please enquire for pricing).

All-inclusive KASP service guidance notes

#### 2. Project workflow

Before commencing their all-inclusive service project, the customer will have a detailed discussion with a Biosearch Technologies representative (e.g. Client Executive, Sales Representative), to ensure all documentation, sample and shipping requirements are met.

Below outlines the details of the processes involved in the KASP service:



#### 3. Sample and sequence submission

Please see the tables in Appendix A, B and C for information on DNA and tissue submission requirements.

It is important to ensure that all BioArk Sample Collection Kits are filled to capacity with 95 samples, leaving one well empty for the Biosearch Technologies non-template control (NTC). Should partially filled kits be submitted where there were sufficient samples to fill completely, extra costs may be chargeable.

#### 4. Data analysis and results

 Following completion of the laboratory work, all results will be returned in Excel format, which can be viewed through our <u>SNPViewer</u> or <u>KlusterCaller</u> software.

All-inclusive KASP service guidance notes

### 5. Retention and return of project samples

- All project material (samples, DNA and data) will be stored for up to 3 months, after the completion of the project. After 3 months, all project materials will be disposed, as per our standard Terms and Conditions.
- Should longer-term storage be required or have any material returned, please discuss this with your Project Manager.
- Additional services outside of the all-inclusive service are available. Please enquire with your Project Manager.

### 6. Useful links and contacts

- Please confirm with your Project Manager the address of where to dispatch your samples.
- For all-inclusive service related queries, please contact <a href="mailto:servicelab-uk@lgcgroup.com">servicelab-uk@lgcgroup.com</a>

All-inclusive KASP service guidance notes

### Appendix A: Forms for services

		Forms to be completed and	I submitted for service project		New Biosearch Tecl	nnologies customer
Application	Biosearch Technologies service	If submitting tissue samples	If submitting DNA	Other useful information	Complete form if you do not have a Biosearch Technologies customer number	Complete form if have not previously registered with our Sequencing Shop
Genotyping (with or without DNA purification)	KASP genotyping or All-inclusive ser- vices	SNP submission template	DNA plate template (96-well plate) DNA plate template (384-well plate) DNA list template (for both 96- and 384-well plate) SNP submission template	DNA requirements factsheet How to prepare and send DNA samples Plant species list KASP anchoring document All-inclusive guidance notes KASP assay design	New account application form	n/a

\*For selected services, project results can only be accessed through our Biosearch Technologies Sequencing Shop. Please ensure you have registered with our Sequencing Shop in order to access your data.

All-inclusive KASP service guidance notes

### Appendix B: DNA submission requirements

		Minimum DNA requirements					Preferred and ac	cepted plastic-v	ware for shipping	Shipping requirements	
Application	Biosearch Technologies service	Minimum sample number per batch	Concentration (per sample)	Volume (per sample)	Diluent requirements	Preferred/ accepted nucleic acid assessment	Plate/tube types	Plate seals	Plate/tube labelling	Packing and dispatch conditions	Other information
Genotyping	KASP genotyping	48 (with 1 well unused for non-tem- plate con- trols) All-inclu- sive: 48	Dry DNA 10 ng per SNP* Wet DNA 7.5 ng per SNP, minimum concentration of 5 ng/µL*	<u>Dry DNA</u> n/a <u>Wet DNA</u> 20 μL**	10 mM Tris- HCI, 0.1 mM EDTA, pH 8.3	Fluorimetric (preferred), spectrophotome tric (accepted)	96-well PCR compatible plates 96-well v-bot- tom plates 96, deep-well 0.5 mL plates (CoStar or Greiner) 384, square- well plates (Greiner, Abgene or Nunc)	Heat-seal- ing (pre- ferred) Cap-mat sealing (ac- cepted)	Barcode (pre- ferred) Name length: 1-15 charac- ters, only al- phanumeric and under- score. If send- ing more than 10 plates, pre- fix with leading zero (e.g. Plate_01) (ac- cepted) If well location identifiers are hidden, please idenitfy well A1	Wet DNA Frozen (wet or dry-ice) Dry DNA Ambient	Plates sizes should not be ≥0.8 mL.

\*Based on human genome-sized organism. Higher DNA concentrations will be required for larger genomes. Please refer to DNA requirements factsheet (see Appendix A) for further guidance. \*\*Volume of DNA required will vary depending on DNA concentration and scope of project (e.g. number of samples/SNPs). Please confirm volumes with your Project Manager before shipping. gDNA = genomic DNA; HMW = high molecular weight

All-inclusive KASP service guidance notes

### Appendix C: Tissue submission requirements

Please note the following important information:

- No biological tissue material (of any type) should be shipped in 0.2 mL tubes or plates
- No GMO seed or BSL-III material is accepted at any Biosearch Technologies site. Non-viable GMO material (e.g. leaves) is accepted.
- No BSL-II material is accepted at the Middleton, WI, USA site. If you wish to process this material, please contact <u>techsupport@lgcgroup.com</u> for alternative options.
- Overall genotyping quality and success rate can be detrimentally affected by poor quality starting material. Biosearch Technologies cannot accept responsibility for lower genotyping call rate/SNP success rate if starting material is not optimal.
- Any deviation from standard sample shipment guidelines (without prior discussion with Biosearch Technologies) may result in refusal of sample receipt or incur additional costs

All-inclusive KASP service guidance notes

### Tissue submission requirements table

		Minimum sample requirements			Preferr	ed and accepted plas	tic-ware for shi	Shipping requirements		
Application	Sample type	Sample number per batch	Volume/ quantity	Quality	Recommended sample collection tube/plate (NOT for submission)	Plate/tube types for submission	Plate seals	Sample handling	Packing and dispatch conditions	Other information
Purification or Purifi- cation and geno-	Leaves	48 or 95 (with at least 1 well unused for non-template controls)*1	4-9, 6 mm (diameter) leaf discs* <sup>2</sup> Leaf mate- rial must be in discs or pre-cut into small pieces, oth- erwise addi- tional pro- cessing costs will be incurred	Young leaf tis- sue, L3/L4 stage, avoiding midrib and axil- lary veins, consistent positions across plants	<u>BioArk Leaf Kit</u> (preferred)	BioArk Leaf Kit (Biosearch Tech- nologies KBS- 9370-001-L) (pre- ferred) 96-well, 2.2 mL plate (preferred for all sites ex- cluding Berlin) 96-well, 1.2 mL plate (preferred for Berlin, ac- cepted for all other sites)	BioArk Leaf Kit (preferred)	If not using the BioArk Leaf Kit, samples should be stored: At -20 °C (KASP genotyping only) or dried at room temperature	BioArk Leaf Kit Ambient Dried Ambient -20 °C Frozen Ethanol Ambient	If HMW DNA is required for other downstream applications (e.g. NGS), it is recommended to immedi- ately dry leaves after harvesting (freezing of fresh leaves will result in DNA degradation)
(all-in- clusive)	Seeds	48 or 95 (with at least 1 well unused for non-template controls)*1	Dependent on size of project - please dis- cuss with your Project Manager	n/a	<u>BioArk Seed Kit</u> ( <u>preferred)</u>	BioArk Seed Kit (Biosearch Tech- nologies KBS- 9370-001-S) (pre- ferred) 96-well, 2.2 mL Porvair plates (preferred) 96-well, 1.2 mL plates (accepted)	BioArk Seed Kit (preferred)	No special sample handling	<u>BioArk Seed Kit</u> Ambient	n/a

All-inclusive KASP service guidance notes

#### Tissue submission requirements table continued

		Minimum	sample requirer	ments	Preferr	ed and accepted plast	tic-ware for shi	Shipping requirements		
Application	Sample type	Sample number per batch	Volume/ quantity	Quality	Recommended sample collection tube/plate	Plate/tube types	Plate seals	Sample handling	Packing and dispatch conditions	Other information
	Blood	48 or 95 (with at least 1 well unused for non-template controls)*1	200 µL	Fresh frozen EDTA, citrate or heparin	EDTA vacutainer PAXgene DNA Blood Tube (BD Bioscience 761165)	96-well, 2.2 mL plate (Porvair Sci- ences 219030) (preferred) 96-well, 1.2 mL plates (accepted)	Heat- sealng (preferred)	No special sample handling	<u>-20 °C</u> Frozen on dry-ce	n/a
Purification or Purifi- cation		48 or 95 (with at least 1 well unused for non-template controls)*1	500 µL	Fresh frozen EDTA, citrate or heparin	EDTA vacutainer PAXgene DNA Blood Tube (BD Bioscience 761165)	96-well, 2.2 mL plate (Porvair Sci- ences 219030) (preferred)	Heat- sealng (preferred)	No special sam- ple handling	-20 °C Frozen on dry-ce	n/a
and geno- typing (all-in- clusive)	Buccal swabs	48 or 95 (with at least 1 well unused for non-template controls)*1	One swab/sample	n/a	Genotek Collec- tion Kit (Oragene OG-500) DNA Buccal Swabs (Isohelix SK-1S)	Swabs Sorted into boxes or bags <u>Tubes</u> Sorted into boxes <u>Discs</u> Sorted into bags	n/a	No special sample handling	Ambient	n/a
	Saliva	48 or 95 (with at least 1 well unused for non-template controls)*1	One swab/One tube or disc/sample	n/a	Genotek Collection Kit (Oragene OG- 500)	Tubes Sorted into boxes Discs Sorted into bags	n/a	No special sample handling	Ambient	n/a

<sup>\*1</sup> If shipping samples to our Berlin, Germany site, a minimum sample number of 24 is required for purification projects.

<sup>12</sup> The exact number of leaf discs will depend on various factors (age of leaf, size/scale of project, plant species etc.). For the exact number of leaf discs for your project, please discuss with your Project Manager. HMW = high molecular weight

The quality management system of LGC UK Ldt. in Hoddesdon, UK is certified by DIN EN ISO 9001:2015.

All orders are subject to LGC, Biosearch Technologies Standard Terms and Conditions.

### Integrated tools. Accelerated science.

**f in** @LGCBiosearch

biosearchtech.com

All trademarks and registered trademarks mentioned herein are the property of their respective owners. All other trademarks and registered trademarks are the property of LGC and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details. No part of this publication may be reproduced or transmitted in any form or by any means, electronic or mechanical, including photocopying, recording or any retrieval system, without the written permission of the copyright holder. © LGC Limited, 2020. All rights reserved. GEN/0356/MW/0520

