

IUB Base Code Guide

Standard Bases

Code	Base
A	Adenosine
C	Cytidine
G	Guanosine
T	Thymidine

Mixed Bases (Wobble)

Code	Base
B	C, G, or T
D	A, G, or T
H	A, C, or T
V	A, C, or G
R	A or G
Y	C or T
K	G or T
M	A or C
S	G or C
W	A or T
N	Any Base



2199 South McDowell Blvd, Petaluma, CA 94954-6904 USA
+1.800.GENOME.1 | +1.415.883.8400 | +1.415.883.8488 (fax)
info@biosearchtech.com

www.biosearchtech.com