

Shipping instructions

For single sample sequencing services

- Transfer your DNA samples into individual 1.5 mL tubes.
- Label each DNA sample with one barcode label
- Close the lid tightly. Please do not use parafilm or adhesive tape to seal the tubes.
- Ship the samples to LGC in a padded envelope (preferably using one of our sample bags)
- If you are supplying your own sample-specific primer, please label the tube with the primer name. LGC will label all submitted primers with a specific barcode upon receipt.

Please send all samples to:

LGC Genomics GmbH
Postfach 940327
12443 Berlin
Germany

Tube type

1.5 mL tube
(Safe-Lock
acceptable)



0.5mL



Screw cap



Different tube




Barcode label

Label horizontal



Label vertical



Tube covered
with tape



Label on
plastic bag



For information on the requirements for our other sequencing services, please see our online sequencing ordering guide (<http://www.lgcgenomics.com/onlineshopmanual/>).

For further information please contact us on

Tel: +49 (0)30 5304 2230 or sequencing@lgcgenomics.com

www.lgcgenomics.com

Germany • United Kingdom • United States